


DEQAS July 2017 - 1,25-dihydroxyvitamin D Method Means * (+/- 1 SD)



	Sample 406	Sample 407	Sample 408	Sample 409	Sample 410
■ DiaSorin Liaison XL (n=98)	140.8	80.3	94.1	117.3	94.2
■ DIASource CT (n=5)	159.8	101.9	95.8	123.2	104.0
■ IDS EIA (n=6)	148.2	89.0	108.3	126.7	102.8
■ IDS RIA (n=15)	181.2	106.1	126.7	141.6	120.0
■ IDS iSYS (n=14)	183.1	106.6	121.3	142.0	121.3
■ IDS iSYS New (n=14)	115.7	73.4	88.5	94.5	79.3
■ LC-MS/MS (n=13)	143.8	93.2	97.0	112.0	99.7
■ ALTM (n=165)	144.7	84.7	98.4	118.9	97.1