


DEQAS July 2018 - 1,25-dihydroxyvitamin D Method Means * (\pm 1 SD)



	Sample 426	Sample 427	Sample 428	Sample 429	Sample 430
■ DiaSorin Liaison XL (n=97)	48.9	103.3	74.3	98.5	97.7
■ DIASource CT (n=5)	62.9	107.0	105.0	115.6	118.1
■ IDS EIA (n=6)	45.5	100.5	99.9	105.6	109.9
■ IDS RIA (n=11)	58.0	113.9	110.5	124.5	124.3
■ IDS iSYS (n=10)	58.8	118.0	123.4	144.5	138.6
■ IDS iSYS New (n=19)	48.7	91.1	92.0	100.7	103.6
■ LC-MS/MS (n=14)	53.5	104.9	98.2	105.1	105.4
■ ALTM (n=162)	50.4	103.3	83.8	103.3	103.3

* Methods with 5 or more results returned