


DEQAS January 2017 - 25-OHD Method Means (+/-1SD) for Major Method Groups*



	Sample 506	Sample 507	Sample 508	Sample 509	Sample 510
■ Abbott Architect New Kit (n=62)	54.3	76.3	28.9	70.8	147.5
■ Abbott Architect Old Kit (n=21)	54.9	74.3	28.6	67.9	146.8
■ Beckman Unicel (n=33)	42.9	63.3	26.4	56.6	141.0
■ DiaSorin Liaison (n=225)	54.2	70.9	29.3	64.9	127.5
■ Euroimmun ELISA (n=11)	48.6	73.8	34.4	65.1	137.0
■ IDS EIA (n=13)	59.8	71.9	34.1	70.0	122.1
■ IDS iSYS (n=80)	54.6	69.9	30.4	70.4	120.9
■ Roche Total (n=170)	44.8	75.0	26.1	65.4	140.5
■ Siemens (n=67)	54.7	74.8	31.5	72.9	114.3
■ HPLC (n=22)	56.1	74.5	31.4	75.9	146.1
■ LC-MS/MS (n=160)	56.2	75.8	31.7	73.1	139.8
■ ALTM (n=915)	52.4	73.1	29.4	67.9	133.6
■ TARGET VALUE	55.8	74.6	31.4	71.6	134.6

* Methods with 10 or more results returned